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- 1) Baker M. Biorepositories: Building better biobanks. *Nature*. 2012;486(7401):141-146. <https://doi.org/10.1038/486141a>.
- 2) Kawasaki H, Itoh T, Takaku Y, *et al.* The NanoSuit method: a novel histological approach for examining paraffin sections in a nondestructive manner by correlative light and electron microscopy. *Lab Invest*. 2020;100(1):161-173. <https://doi.org/10.1038/s41374-019-0309-7>.
- 3) Takaku Y, Suzuki H, Ohta I, *et al.* A thin polymer membrane, nano-suit, enhancing survival across the continuum between air and high vacuum. *Proc Natl Acad Sci U S A*. 2013;110(19):7631-7635. <https://doi.org/10.1073/pnas.1221341110>.

- 4) Takaku Y, Suzuki H, Kawasaki H, *et al.* A modified 'NanoSuit(R)' preserves wet samples in high vacuum: direct observations on cells and tissues in field-emission scanning electron microscopy. *R Soc Open Sci.* 2017;4(3):160887. <https://doi.org/10.1098/rsos.160887>.
- 5) Ilgen P, Stoldt S, Conradi LC, *et al.* STED super-resolution microscopy of clinical paraffin-embedded human rectal cancer tissue. *PLoS One.* 2014;9(7):e101563. <https://doi.org/10.1371/journal.pone.0101563>.
- 6) Creech MK, Wang J, Nan X, Gibbs SL. Superresolution Imaging of Clinical Formalin Fixed Paraffin Embedded Breast Cancer with Single Molecule Localization Microscopy. *Sci Rep.* 2017;7:40766. <https://doi.org/10.1038/srep40766>.
- 7) Shinmura K, Kusafuka K, Kawasaki H, *et al.* Identification and characterization of primary cilia-positive salivary gland tumours exhibiting basaloid/myoepithelial differentiation. *J Pathol.* 2021;254(5):519-530. <https://doi.org/10.1002/path.5688>.
- 8) Shinmura K, Kato H, Kawasaki H, *et al.* Primary Cilia Are Frequently Present in Small Cell Lung Carcinomas but Not in Non-Small Cell Lung Carcinomas or Lung Carcinoids. *Lab Invest.* 2023;103(2):100007. <https://doi.org/10.1016/j.labinv.2022.100007>.
- 9) Itoh T, Yamada S, Ohta I, *et al.* Identifying Active Progeny Virus Particles in Formalin-Fixed, Paraffin-Embedded Sections Using Correlative Light and Scanning Electron Microscopy. *Lab Invest.* 2023;103(1):100020. <https://doi.org/10.1016/j.labinv.2022.100020>.
- 10) Yamada S, Itoh T, Ikegami T, *et al.* Association between human papillomavirus particle production and the severity of recurrent respiratory papillomatosis. *Sci Rep.* 2023;13(1):5514. <https://doi.org/10.1038/s41598-023-32486-8>.
- 11) Scimeca M, Pietroiusti A, Milano F, *et al.* Elemental analysis of histological specimens: a method to unmask nano asbestos fibers. *Eur J Histochem.* 2016;60(1):2573. <https://doi.org/10.4081/ejh.2016.2573>.
- 12) Shinmura K, Kawasaki H, Baba S, *et al.* Utility of Scanning Electron Microscopy Elemental Analysis Using the 'NanoSuit' Correlative Light and Electron Microscopy Method in the Diagnosis of Lanthanum Phosphate Deposition in the Esophagogastrroduodenal Mucosa. *Diagnostics (Basel).* 2019;10(1). <https://doi.org/10.3390/diagnostics10010001>.
- 13) Ooishi M, Yamada S, Itoh T, *et al.* Diagnosis of Ion-Exchange Resin Depositions in Paraffin Sections Using Corrective Light and Electron Microscopy-NanoSuit Method. *Diagnostics (Basel).* 2021;11(7). <https://doi.org/10.3390/diagnostics11071193>.
- 14) Sugiura K, Kawasaki H, Egami T, *et al.* Iron deposition in gastric black spots: Clinicopathological insights and NanoSuit-correlative light and electron microscopy analysis. *DEN Open.* 2024;5(1):e398. <https://doi.org/10.1002/deo2.398>.